## GreenGel DNA/RNA Stain

### Datenblatt / Data sheet

<b>Art.No. BS30.200.0050</b>	<b>50 µl</b>	<b>(ca. 15 Mini Gele)</b>
<b>Art.No. BS30.200.0500</b>	<b>500 µl</b>	<b>(ca. 150 Mini Gele)</b>
<b>Art.No. BS30.200.1000</b>	<b>1000 µl</b>	<b>(ca. 300 Mini Gele)</b>

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## Description

The Bio&SELL GreenGel DNA/RNA Stain is suitable for the (simultaneous!) staining of DNA (double- and single-stranded) as well as RNA and can be used in agarose gels as well as in polyacrylamide gels. The sensitivity is thereby **up to 4 x higher** compared to ethidium bromide, so that amounts from approx. 0.1 ng DNA can be detected. The staining solution is added directly into the gel and into the running buffer, destaining after electrophoresis is omitted; this saves valuable time. The stimulation of the Bio&SELL GreenGel DNA/RNA Stain is performed with blue light (wavelength 470 - 490 nm), detection with SYBR® Green filters. A special feature is the two emission peaks of the dye: DNA bands appear green, RNA bands red.


Fig. : Absorption and emission spectrum of GreenGel DNA/RNA Stain for DNA.  
Excitation maximum: 490 nm (blue light)  
Emission maxima: 520 nm for DNA, 635 nm for RNA (not shown)

Fig. Top right: Agarose gel stained with GreenGel DNA/RNA Stain

**The Bio&SELL GreenGel DNA/RNA Stain is thus a more comfortable, more sensitive, safer and environmentally friendly alternative to the conventional ethidium bromide.**

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--

## Applications:

Staining of DNA (double- and single-stranded) as well as RNA in agarose and polyacrylamide gels

## Reaction Conditions:

- 1x solution directly in the gel
- 0.5x solution in the electrophoresis buffer
- 5 µl GreenGel DNA/RNA Stain per 100 ml of agarose or polyacrylamide gel
- 2.5 - 3 µl GreenGel DNA/RNA Stain per 100 ml 1x electrophoresis buffer

## Note:

- Destaining of the gel is not required, but in exceptional cases it can improve the visibility of the band by reducing the background.
- Staining of the gel after the run is not recommended.
- Do not expose the gels to ultraviolet light for too long (especially for cloning).
- Replace the electrophoresis buffer after a few runs.
- If the gel solution is boiled several times, add 2 - 3 µl of GreenGel DNA/RNA Stain per boil.

## Protocol for agarose-geles:

1. Boil the agarose solution as usual.
2. Let the solution cool down until a temperature of 50 -60 ° C is reached.
3. Immediately before pouring the gel, add 5 µl GreenGel DNA/RNA Stain per 100 ml gel solution and mix it evenly, but gently, so that no air bubbles are formed.
4. Pour the gel as usual and use the combs.
5. Prepare the required volume of 1x running buffer.
6. Add 2.5 - 3 µl GreenGel DNA/RNA Stain per 100 ml electrophoresis buffer.
7. Place running buffer in the chamber and perform electrophoresis as usual.
8. The stimulation is performed with blue light (wavelength 470 - 490 nm), detection with SYBR® Green filters

## Protocol for polyacrylamid-geles:

1. Prepare the polyacrylamid solution as usual.
2. Add TEMED and APS and immediately add **5 µl GreenGel DNA/RNA Stain per 100 ml gel.**
3. Mix it evenly, but gently, so that no air bubbles are formed.
4. Pour gel as usual and use combs.
5. Prepare the required volume of 1 x running buffer
6. **Add 2.5 – 3 µl GreenGel DNA/RNA Stain per 100 ml electrophoresis buffer.**
7. Add running buffer to the chamber, perform electrophoresis as usual.
8. The stimulation is performed with blue light (wavelength 470 - 490 nm), detection with SYBR® Green filters

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--

# GreenGel DNA/RNA Stain

## Technical Data

**Shipment:** at room temperature

**Storage:** **light-protected at +4 ° Celsius for 24 months**

### Safety and disposal:

This product should only be used by persons who have routine in laboratory applications. Wear appropriate protective clothing such as gowns, gloves and goggles. In case of contact with skin and eyes the affected areas should be immediately washed with water or rinsed out.

The non-carcinogenicity of GreenGel DNA/RNA Stain was detected in the Ames test. Microcore tests on erythrocytes from mouse bone marrow and chromosome aberration tests on mouse spermatocytes gave no evidence of mutagenicity.

Despite their proven lower mutagenicity, all nucleic acid dyes should be handled with care. Used running buffers containing dye should be cleaned as well as molten gels using a suitable filter (e.g., activated carbon). As soon as fluorescence is no longer detectable, the liquids can be disposed of with plenty of water in the drain. Alternative is the autoclaving of all used fluids since the GreenGel DNA/RNA Stain degrades during prolonged exposure at temperatures above 100 ° C.

0.5x concentrated solutions of the GreenGel DNA/RNA Stain (electrophoresis buffer) can be disposed of with sufficient water in the effluent.

Used gels can be disposed of with the laboratory waste.

Contact your safety officer for current local regulations, as they may differ or are changed from this recommendation.

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--

## Bio&SELL Agarose products you may also be interested in:

### Bio&SELL Universal Agarose

Bio&SELL Universal Agarose is a cost-saving standard agarose ideal for analytical and preparative separation of nucleic acids in the range of 0.05 - 50 kbp.

€ 189,-/ 500g

### Bio&SELL Low Melt Agarose

Bio&SELL Low Melt Agarose is a high-quality agarose with a relatively low melting point, ideal for the separation of nucleic acids in the range of 0.8 – 25 kbp. Due to its high-resolution capabilities it is suitable for all preparative gel electrophoresis.

€ 94.90/ 25g

### Bio&SELL High resolution Agarose

Bio&SELL High-definition Agarose as an alternative for the work intensive, expensive polyacrylamid gels. It is optimized for the separation of nucleic acids in the range of 20 bp – 800 bp.

€ 89,-/ 25g

**The agaroses with the best price performance ratio!**

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--

## GreenGel DNA/RNA Stain

**Bio&SELL DNA Ladders you may also be interested in:**

**Sharp bands and excellent resolution in every size range!**


**Choose your Marker:**

For the size range 50 bp to 700 bp: Bio&SELL 50 bp DNA marker

For the size range 100 bp to 1500 bp: Bio&SELL 100 bp DNA marker

For the size range 250 bp to 10 000 bp: Bio&SELL 1 kb DNA marker

For the size range 125 bp to 23 130 bp: Bio&SELL  $\lambda$ -HindIII DNA marker

**Available with or without loading dye!**

## Products for cDNA Synthesis you may also be interested in: Reverse Transcriptases

- **Simple handling of the kits**
- **High reproducibility**
- **Fast protocols**

**For high quality cDNA perfectly fitted for further applications.**

### **SCRIPTUM-37 reverse transcriptase**

Particularly suitable for routine applications with simple templates.

Temperature optimum of this reverse transcriptase is 37 °C.

99.90 €/ 10,000 units

### **SCRIPTUM-55 reverse transcriptase**

For simple and difficult templates.

Optimum temperature for the synthesis of longer cDNA fragments (>7kb) is 50 °C.

139.90 €/ 10,000 units

### **SCRIPTUM-55 reverse transcriptase „improved“**

Specifically designed for difficult templates.

Increased thermal stability: working temperature up to 55 °C

Highly sensitive and highly processive!

249.90 €/ 10,000 units

**Bio&SELL reverse transcriptases:**

**Specific  
High cDNA yields  
Effective**

**[www.bio-sell.com](http://www.bio-sell.com)**

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--