HANDBUCH / MANUAL

QPCR Mixe EvaGreen®

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Telefon : +49 (0) 9128 – 724 32 32
E-Mail: info@bio-sell.de

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Table of content

Protocols for 5x QPCR Mixes EvaGreen® Rox, No Rox, Capillary	p. 15
Protocols for 5x QPCR Mixes EvaGreen® ROBUST	p. 16
Additional protocol for optional UNG treatment only for ROBUST mix	p. 17
Cycler compatibility	p.18

The 5x QPCR Mixes EvaGreen® at a glance

The 5x QPCR Mixes EvaGreen® are optimized "ready-to-use" solutions for real-time applications such as quantitative PCR (QPCR), which require the EvaGreen® or SYBR® Green dye. They contain all the necessary components for a quantitative PCR except template DNA and primers.

EvaGreen®- dye

The modern dye EvaGreen® is compatible with every thermal cycler and real-time PCR protocol, simply choose the default setting for SYBR® Green or FAM.

EvaGreen® – Advantages at a glance

- Lower PCR inhibition than SYBR® Green - can be used in higher concentrations
- Higher sensitivity - Signal stronger than SYBR® Green
- Unbeatable high-resolution DNA melting curves
- Thermal stability – comfortable handling
- Non-cytotoxic, non-mutagenic
- Suitable for all QPCR devices

Choose between four different QPCR EvaGreen® Mixes:

Rox-Mix: is compatible with all real-time cyclers that need Rox as an internal reference. (Rox. passive internal reference dye, which is used for normalization of fluorescence reporter signal during the QPCR).

No Rox-Mix: is compatible with real-time cyclers that do not require a reference dye.

Capillary-Mix: is specifically designed for real-time cycler who work with capillaries and do not require a reference dye.

ROBUST-Mix: the "ROBUST" Mix is, due to the specific buffer composition and the GC enhancers, in particular suitable for use in GC-rich templates and templates with long and frequent repeats. Also recommended for UNG treatment (protocol p.17). The 5x EvaGreen® QPCR "ROBUST" is designed for use in 'standard cycling mode' on 'Standard' and 'Fast' real-time cyclers, regardless of their need of a reference dye.

Summary of the cycler and probe compatibility of 5x QPCR Mixes EvaGreen® on p.18.

Application examples of the 5x QPCR Mixes EvaGreen®

- Detection and quantification of DNA and cDNA targets
- Profiling gene expression
- Microbial detection
- Determination of viral loads
- High-resolution melting curve analysis (HRM)
- UNG treatment (here we recommend the 5x QPCR Mix EvaGreen® (ROBUST))

Composition

- Hot-Start Taq-DNA-polymerase
- 5x EvaGreen® QPCR reaction buffer
- 12.5 mM MgCl₂, (= 2.5 mM in 1x solution)
- dATP, dGTP, dCTP and dTTP in an optimized composition
- EvaGreen®-dye
- BSA (only included in the 5x QPCR Mix EvaGreen® Capillary)
- Rox-dye (not included in the 5x QPCR Mixes EvaGreen® „No Rox“ and „Capillary)

Additionally included in the 5x QPCR Mix EvaGreen® (ROBUST):

- dNTP mix with dUTP/ dTTP
- GC enhancer
- Blue dye for visualization

5x QPCR Mixes EvaGreen® (Rox, No Rox, Capillary)

Recommended pipetting protocol

The single components should be pipetted as follows (here for 20 µl end volume):

Component	Volume	Final concentration
5x QPCR Mix EvaGreen®	4 µl	1 x
Primer 1	x µl	80 - 250 nM
Primer 2	x µl	80 - 250 nM
Template-DNA	1 - 5 µl	1 -50 ng/ µl
Water	Fill up to 20 µl	

Recommended PCR protocol

To activate the polymerase in the 5x QPCR mix EvaGreen®, include an incubation step for 15 minutes at 95°C at the beginning of the QPCR cycle.

Cycle step	Temperature	Time	Number of cycles
Initial Incubation	95 °C	15 min	1 x
Denaturation	95 °C	15 sec	40 x
Annealing	60 - 65 °C	20 sec	40 x
Elongation	72 °C	20 sec	40 x

5x QPCR Mixes EvaGreen® (ROBUST)

Recommended pipetting protocol

The single components should be pipetted as follows (here for 20 µl end volume):

Component	Volume	Final concentration
5x QPCR Mix EvaGreen® (ROBUST)	4 µl	1 x
Primer 1 (10 pmol/ µl)	0.2 - 0.4 µl	100 - 200 nM
Primer 2 (10 pmol/ µl)	0.2 - 0.4 µl	100 - 200 nM
Template-DNA	1 - 5 µl	0.002 - 2 ng/ µl
Water	Fill up to 20 µl	

Recommended PCR protocol

To activate the polymerase in the 5x QPCR mix EvaGreen®, include an incubation step for 12 minutes at 95°C at the beginning of the QPCR cycle.

Cycle step	Temperature	Time	Number of cycles
Initial Incubation	95 °C	12 min	1 x
Denaturation	95 °C	15 sec	40 x
Annealing	60 - 65 °C	20 - 30 sec	40 x
Elongation	72 °C	20 - 30 sec	40 x

Additional protocol for optional UNG treatment only for 5x QPCR Mix EvaGreen® (ROBUST)

Recommended pipetting protocol

The single components should be pipetted as follows (here for 20 µl end volume):

Component	Volume	Final concentration
5x QPCR Mix EvaGreen® (ROBUST)	4 µl	1 x
Primer 1 (10 pmol/ µl)	0.2 - 0.4 µl	100 - 200 nM
Primer 2 (10 pmol/ µl)	0.2 - 0.4 µl	100 - 200 nM
UNG (Uracil-N-Glykosylase)	x µl *	0.01 U/ µl
Template-DNA	1 - 5 µl	0.002 - 2 ng/ µl
Water	Fill up to 20 µl	Template-DNA

* according to the manufacturer's instructions

Recommended PCR protocol

To activate the polymerase in the 5x QPCR Mix EvaGreen®, include an incubation step for 12 minutes at 95°C at the beginning of the QPCR cycle.

Cycle step	Temperature	Time	Number of cycles
UNG treatment	50 °C	2 min	1 x
Initiale Incubation	95 °C	12 min	1 x
Denaturation	95 °C	15 sec	40 x
Annealing	60 - 65 °C	20 - 30 sec	40 x
Elongation	72 °C	20 - 30 sec	40 x

Cycler compatibility:

QPCR Mix Rox	QPCR Mix No Rox	QPCR Mix Capillary	QPCR Mix ROBUST
ABI PRISM® 5700, 7000,7300, 7500, 7700, 7900, 7900HT, (inkl. Fast- Block)	Corbett Rotor-Gene™ 3000	Roche Light Cycler® 1.x	Applied BioSystems: QuantStudio™ 12K Flex, ViiA™ 7, 7900HT,ABI Prism® 7500,7700 SDS, StepOne™, StepOne™ Plus
StepOne™, Step One™ Plus	Corbett Rotor-Gene™ 6000	Roche LightCycler® 2.	Bio-Rad: CFX96™, CFX384™ iQ™5, MyiQ™, Chromo™ 4, Opticon™ 2, MiniOpticon™
Stratagene: Mx3000P™, Mx3005P™, Mx4000®	Eppendorf Mastercycler® ep realplex		Bio-Rad: CFX96™, CFX384™, iQ™5, MyiQ™, Chromo™ 4, Opticon™2, MiniOpticon™
	Roche Light- Cycler® 480		Corbett: Rotor-Gene™ Q, Rotor-Gene™ 6000
	Bio-Rad/ MJ Research: Opticon™, Opticon™ 2, Mini Opticon, CFX 96/ 384, Chromo 4™, iQ5™		Eppendorf: Mastercycler® ep, realplex2, ep realplex4
			Illumina: Eco™
			Roche: LightCycler®480

Technical Data:

Shipment: at room temperature

Storage: -20 °Celsius

Storage at room temperature for up to 1 month has no detrimental effects on the quality of the QPCR mix EvaGreen®.

Safety warning: This product and its components should be handled only by persons trained in laboratory techniques. It is advisable to wear suitable protective clothing, such as laboratory overalls, gloves and safety glasses. Care should be taken to avoid contact with skin or eyes. In case of contact with skin or eyes, wash immediately with water.

Note

Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required

EvaGreen® is a registered trademark of BIOTIUM, Inc., Hayward, CA, USA.

SYBR® is a registered trademark of Molecular Probes, Inc.

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer, where such research does not include testing, analysis or screening services for any third party in return for compensation on a per test basis. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes.

Product	Content	Order No.	Price
5x QPCR Mix EvaGreen® (Rox)	200 µl = 50 reactions	BS 76.580.0200	24.90 €
	1000 µl = 250 reactions	BS 76.580.1000	119 €
	5000 µl = 1250 reactions	BS 76.580.5000	535 €
5x QPCR Mix EvaGreen® (No Rox)	200 µl = 50 reactions	BS 76.590.0200	24.90 €
	1000 µl = 250 reactions	BS 76.590.1000	119 €
	5000 µl = 1250 reactions	BS 76.590.5000	535 €
5x QPCR Mix EvaGreen® (Capillary)	200 µl = 50 reactions	BS 76.570.0200	24.90 €
	1000 µl = 250 reactions	BS 76.570.1000	119 €
	5000 µl = 1250 reactions	BS 76.570.5000	535 €
5x QPCR Mix EvaGreen® (ROBUST)	200 µl = 50 reactions	BS 76.550.0200	39 €
	1000 µl = 250 reactions	BS 76.550.1000	174 €
	5000 µl = 1250 reactions	BS 76.550.5000	779 €

You are looking for information on the quantitative real-time-PCR (QPCR) or plan to establish the method in your laboratory?

Bio&SELL has the right product for you!

- ✓ 5x QPCR Mixes („ready-to-use“)
- ✓ Labeled probes (taqman® probes, molecular beacons®)
- ✓ Seminars about real-time-PCR
- ✓ Real-time-PCR plastic
- ✓ One-Step RevTrans-QPCR Kits
- ✓ RNA isolation kits

Find out more about the wide product portfolio of QPCR products of Bio&SELL on **www.bio-sell.com** or contact us directly!

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