# Freezeprotect

## Datenblatt / Data sheet

<b>Art.No. BS33.200.050</b>	<b>50 ml</b>
<b>Art.No. BS30.200.100</b>	<b>100 ml</b>
<b>Art.No. BS30.200.500</b>	<b>500 ml</b>

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## Freezprotect

### Description:

The serum-free, ready-to-use freezing medium **Freezprotect** from Bio&SELL is optimized for freezing animal and human cells (adherent and suspension cells). It contains no human or animal components and is especially recommended for cells from serum-free culture.

### Application and special advantages:

- For cryopreservation of human and animal cells
- Due to its serum-free formulation, it is particularly suitable for the preservation of **serum-free cultured cells**.
- The specially optimized composition of the medium guarantees a **very high cell viability** after the thawing process.
- Suitable **for all cell types**, including primary human cells and human cell lines.
- The risk of contamination with viruses, BSE or interfering proteins is excluded as it contains no human or animal components.

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--

## Freezprotect

### Composition of Freezprotect Freezing Medium:

- Chemically defined and optimized mixture of salts, sugar, DMSO and other antifreeze-substances.
- It contains no animal and human components.

### Technical Data

**Shipment:** room temperature

**Storage:** -20°C (stable for 24 month)

Bio&SELL Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
--	--	--

## Freezing protocol:

**Note:** In general, use only vital cells in the logarithmic growth phase for getting optimal results.

- Thaw the Freezprotect freezing-media and store till use at 2-8 ° C (fridge).
- Trypsinize adherent cells as usual and transfer the cells into cell culture medium. Stop the trypsin activity by adding a trypsin inhibitor.
- Centrifuge the cells as usual.
- Discard supernatant.
- Wash cell pellet in PBS (without Ca<sup>2+</sup> / Mg<sup>2+</sup>).
- Centrifuge again (100 - 200 g, 5 - 10 minutes) and then resuspend the cell pellet in PBS.
- Determine cell number and cell vitality by trypan blue staining.
- Centrifuge again and discard the supernatant.
- Transfer the cells into the cold Freezprotect freezing medium (5x10<sup>5</sup> - 2x10<sup>6</sup> cells / ml Freezprotect).
- Mix the suspension with the cells thoroughly by careful repeated pipetting until no cell clumps are left.

## Freezeprotect

- Transfer the cell suspension into labeled cryotubes (0.5 - 1.5 ml / tube).
- Freeze the cells in an automatically or manually controlled cryo freezer. The optimal freezing rate is about 1 ° C / minute.
- **Alternative:**
  - Place the tubes in the refrigerator for 15 minutes (so that the freezeprotect freezing medium can penetrate into the cells).
  - Then freeze the cryotubes for 2 h at -20 ° C.
  - Add the cryotubes to the vapor phase of liquid nitrogen overnight.
- Store the cryotubes in a cryotank with liquid nitrogen.

## Thaw protocol:

- Remove the cryotubes from the cryotank and thaw them in warm water as fast as possible (<1 minute).
- Disinfect the exterior of the cryovials with alcohol and transfer the cells under sterile conditions to a centrifuge tube containing 10 ml of growth medium and mix gently.
- Centrifuge the cells (150 - 200 g, 5 - 10 minutes).
- Discard the supernatant and recover the cells into the designated culture medium.
- Determine the viability of the cells using a suitable method, such as FACS or Trypan Blue staining.

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