



# High resolution agarose

## Data sheet

<b>Order-No. BS 20.48.0TM</b>	<b>5g</b>
<b>Order-No. BS 20.48.025</b>	<b>25g</b>
<b>Order-No. BS 20.48.100</b>	<b>100g</b>
<b>Order-No. BS 20.48.250</b>	<b>250g</b>

(For research and *in vitro* application only)

---

Lot-No.:

Best before:

Appearance:           homogeneous powder

Colour:                   white

---

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--



# High resolution agarose

## Description

Bio&SELL high resolution agarose is specially designed and optimized for separation of nucleic acids in a range of 20bp – 800bp.

## Analytical specification

Gel strength (1,0%):	> 300 g/cm <sup>2</sup>
Gelling temperature (1,5 %):	< 35 °C
Elektroendosmosis (-mr)	< 0,5
Amount of water	< 10,0 %
DNA-Binding	not detectable
DNase- a. RNase activity	not detectable

## Proposed agarose concentrations

Fragmentsize (bp)	Endconcentration Agarose (%)	
	1x TAE-Puffer	1x TBE-Puffer
50 – 800	2,0	1,8
100 – 600	3,0	2,0
50 – 250	4,0	3,0
20 – 130	5,0	4,0
< 80	-	5,0

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--

# High resolution agarose

## Preparation of the agarose:

▶▶▶ **Note:** For 100 ml of a 1% agarose gel prepare 1 g of agarose in 100 ml electrophoresis buffer.

### Method 1: Microwave

1. Select a flask which can uptake the two- to fourfold of the final volume of gel solution
2. Add **cooled** 1x or 0,5x electrophoresis buffer and a magnetic stir bar to the flask.
3. Place the flask on a magnetic stirrer and slowly intersperse the agarose. Stir fast at all time to avoid clumping.
4. Take out the magnetic stir bar if it is not teflon coated.
5. Incubate the agarose solution for 15 minutes before heating. This will reduce foaming during heating.
6. Weigh the flask.
7. Close the flask with plastic wrap. Make a small opening in the plastic wrap for ventilation.
8. For **agarose concentrations > 4%**, perform the following **additional steps** to avoid foaming:

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--

## High resolution agarose

- A. Heat the flask in the microwave at **medium** power for 1 minute.
  - B. Take out the flask from the microwave.
  - C. Allow the solution to cool down for 15 minutes on the bench.
9. Heat the flask in the microwave at **medium** power for 2 minutes.
  10. Take out the flask from the microwave.  
**Attention: boiling retardation possible!**
  11. **Carfully** shake the flask to resuspend agarose which has not been dissolved yet.
  12. Heat the flask at **high** power in the microwave until the solution is boiling.
  13. **Let boil for 1 minute** until all particles are dissolved.
  14. Take out the flask from the microwave.
  15. **Carfully** shake the flask to mix the solution well.
  16. Weigh the flask again and adjust lost weight by the addition of warm, distilled water.
  17. Mix the solution well.
  18. Let the solution cool down at room temperature until a temperature of 50 – 60 °C is reached. After casting of the gel, let it further cool down at room temperature until the gel is stable.  
**Cool another 20 minutes at 4 °C** to reach optimum resolution and handling.

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--

## High resolution agarose

### Method 2: Heating plate

1. Select a flask which can uptake the two- to fourfold of the final volume of gel solution.
2. Add **cooled** 1x or 0,5x electrophoresis buffer and a magnetic stir bar to the flask.
3. Place the flask on a magnetic stirrer and slowly intersperse the agarose. Stir fast at all time to avoid clumping.
4. Weigh the flask.
5. Close the flask with plastic wrap. Make a small opening in the plastic wrap for ventilation.
6. Allow the solution to boil while stirring it.
7. Allow to boil it further until the agarose is completely dissolved (approx. 10 minutes).
8. Weigh the flask again and adjust lost weight by the addition of warm, distilled water.
9. Mix the solution well.
10. Let the solution cool down at room temperature until a temperature of 50 – 60 °C is reached. After casting of the gel, let it further cool down at room temperature until the gel is stable.  
**Cool another 20 minutes at 4 °C** to reach optimum resolution and handling.

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--



# High resolution agarose

## Technische Daten

**Gel strength of agarose:** Die gel strength, specified in g/cm<sup>2</sup> at a given agarose concentration indicates the stability of the completed gel.

**Gelling temperature of agarose:** gelation range of temperature in which the gel reach its end consistency.

**Elektroendosmosis of agarose:** The smaller the elektroendosmosis the better the agarose. A small value for elektroendosmosis (< 0,12) characterize a high quality agarose, which is most suitable for gel electrophoretic labour.

**Content of sulphate of agarose:** A small content of sulphate indicates a high quality agarose.

---

Shipment: at room temperature

Storage: at room temperature

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--



# High resolution agarose

## Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

<b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg	E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a>	Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33
---	--	--

## High resolution agarose

**Further products which will interest you:**

### Bio&SELL products for cDNA synthesis

**The Bio&SELL SCRIPTUM One-step kits: obtaining high-quality cDNA fast and effective**

**Bio&SELL  
SCRIPTUM Standard Kit**


Combines the specific reverse transcription with the subsequent PCR fast and easy in only one tube.

**Bio&SELL  
SCRIPTUM High Precise**

cDNA synthesis extreme fast and precise because of the Taq-DNA-Polymerase with proof-reading-activity

**Bio&SELL  
SCRIPTUM First (strand)**

For the synthesis of first-strand cDNA: by highly structured and long cDNA fragments, an extreme sensitive and specific RT-PCR and DNA labeling.



**Bio&SELL  
cDNA-Kits SCRIPTUM  
“ready to use”:**

- fast
- effective
- low priced

[www.bio-sell.de/molekularbiologie/reverse-transkription.html](http://www.bio-sell.de/molekularbiologie/reverse-transkription.html)

<p><b>Bio&amp;SELL</b> Lohweg 27 90537 Feucht bei Nürnberg</p>	<p>E-Mail: <a href="mailto:info@bio-sell.de">info@bio-sell.de</a> Internet: <a href="http://www.bio-sell.de">www.bio-sell.de</a></p>	<p>Fon : +49 (0) 9128 – 724 32 32 Fax : +49 (0) 9128 – 724 32 33</p>
--	--	--